

# 3 Simple Steps to Success



**IDEAL PATIENT CRITERIA FOR PRESBYOPES:** Symptomatic presbyopes, current contact lens wearers or those motivated to wear contact lenses, normal cornea and pupil size, astigmatism  $\leq -1.00D$

**1•2•3 See Naturally®** - Follow these best practices for NaturalVue® (etafilcon A) Multifocal 1 Day Contact Lenses when fitting patients.



## 1. STARTING Power



Refract to  
20/15



Use the best corrected  
spectacle refraction  
(BCSR) using the full  
cylinder component



Use the Red/Green  
(Duochrome) test  
for a binocular check  
that the patient is not  
overplussed (one click  
into the Green for best  
starting point)



Then enter the full  
spectacle refraction  
into the NaturalVue®  
Multifocal QuickStart  
Calculator<sup>1</sup>



## 2. SETTLE for 10 minutes

Let lenses settle 10 minutes & distance vision will begin to clear



## 3. Evaluate & Enhance - Distance 1st; then Near

Make an On-Eye Diagnostic Lens Change – The 0.25D Difference to Success:

- **To enhance DISTANCE:** CHANGE THE DIAGNOSTIC LENS POWER on eye by -0.25D (up to -0.50D) in the DOMINANT eye, and -0.25D in the NON-DOMINANT eye only if needed.
- **To enhance NEAR:** CHANGE THE DIAGNOSTIC LENS POWER on eye by +0.25D (up to +0.50D) in the NON-DOMINANT eye, and +0.25D in the DOMINANT eye only if needed.

**For additional fitting support, contact Consultation:  
1-844-VTIVISION (1-844-884-5367), ext. 102**

<sup>1</sup>The NaturalVue Multifocal QuickStart Calculator is available at [www.naturalvuecalculator.com](http://www.naturalvuecalculator.com), or as an app from the Apple® App Store® and Google Play™. Google Play and the Google Play logo are trademarks of Google Inc. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

# NaturalVue® (etafilcon A) Multifocal 1 Day Contact Lenses Vision Enhancement Guide



## For any questions or issues regarding fitting:

- ☑ Check that lenses have settled 10 min
  - ☑ Verify best corrected spectacle refraction used
  - ☑ Verify vertex correction (Round toward minus)
  - ☑ Use Red/Green (Duochrome) test to verify Rx is not overplussed
- For presbyopes, a 0.25D lens change can greatly impact vision. Listen to patient. If the patient says vision is not clear, change the on eye lens power. Make a 0.25D change in the power, one eye at a time, starting with -0.25D in dominant eye to refine distance; +0.25D in non-dominant eye to refine near.
  - Patients in monovision or center near multifocal designs may not have been using their central vision in both eyes for distance. It may take them a few days to fully adapt to this change.

Issue	
Distance vision not clear	<b>Dominant:</b> Change diagnostic lens -0.25 (-0.50 if needed) <b>Non-Dominant:</b> Change diagnostic lens -0.25 only if needed
Distance vision not clear and over refraction does not help	<b>Dominant:</b> Change diagnostic lens -0.25 (-0.50 if needed) <b>Non-Dominant:</b> Change diagnostic lens -0.25 only if needed
Near vision not clear	<b>Non-Dominant:</b> Change diagnostic lens +0.25 (+0.50 if needed) <b>Dominant:</b> +0.25 only if needed
Glare, Halos, vision seems "different"	<b>Dominant:</b> Change diagnostic lens -0.25 (-0.50 if needed) If Center Near MF wearer or monovision, advise may take 2-3 days to adapt
Issues with driving at night	<b>Dominant:</b> Change diagnostic lens -0.25 (-0.50 if needed)
Vision is clear but feels strange, 3D effect, things jumping out at them	<b>Dominant:</b> Change diagnostic lens -0.25 (-0.50 if needed) If Center Near MF wearer or monovision, advise may take 2-3 days to adapt
Have correct Best Corrected Spherical Equivalent and Vertex correction, been adding minus and distance still not clear	<b>Check lenses settled 10 minutes</b> If Center Near MF wearer or monovision, advise may take 2-3 days to adapt <b>Dominant:</b> Change diagnostic lens -0.25 (-0.50 if needed)
Vision fluctuating, going in and out, feels like have to stare at something to focus – Signs of over-minus	<b>Non-Dominant:</b> Change diagnostic lens +0.25 (+0.50 if needed) <b>Dominant:</b> +0.25 only if needed